


Challenges in Transfer of Complex Assay Formats

Cheryl McAlpine, Ph.D.

Merck Millipore, Discovery & Development Solutions
Oxford, UK

Overview

- Assay transfer
 - What & Why?
- Common Challenges
- Best Practice to Approach
- Example Transfer for Cell-based assay
 - Immunogenicity
- Summary & Critical Aspects



Assay Transfer: What & Why?

Transfer of Analytical Method

Starting Point Variable

From partially developed non-validated methodology to fully validated robust procedure

Common Objectives

Demonstration of 'acceptable' assay performance at transfer site

Reasons for Transfer

Transfer from R&D to Regulatory environment

Resource/ technology availability

Leverage cost & efficiency benefits

Accelerate study timelines

Partner with technical expertise

Common Challenges to Assay Transfers

- Design of effective transfer strategy
- Partnership & communication
- Differences between research sites:
 - Equipment
 - Readout technology
 - Reagents
 - Techniques
 - Interpretation of methodology and/or data
- Aggressive timelines

Approach to Assay Transfer

Pre-Transfer

- Early planning
- Understand needs = establish purpose
- Suitability of methods

Develop Effective Documentation

- Transfer study plan: *'What, How, When'*
 - Agreement between sites
 - Detailed, tailored objectives
 - Pre-determined acceptance criteria

Lab Phase

- Method optimised as required for purpose
- Discussion and agreement
- In-line with study plan & recommendations

Reporting

- Final Report of complete study
- Capture solution based approaches
- Demonstrate 'fit for purpose'

Regular Communication

Monitor Constantly

Feedback

Example Transfer

Challenges & Identified Solutions

Cell-Based Assays

- Used in all stages of drug development
- Unique complexity & variability
- Considered most relevant measure of biological mechanism of action



Neutralising Antibody Assay

- Immunogenicity testing

Assay Specific Transfer Challenges

Neutralising Antibody (NAb) Assay

▪ **Variability**

- Target expression levels - adequate, stable, specificity
 - Cell growth conditions/ seeding density
 - Cell passage number

▪ **Limited Dynamic Range**

▪ **Matrix Interference**

- Components may enhance/inhibit activity of drug product in assay

▪ **Choice of signal analysis**

- Understanding **Mechanism of Action**
 - Overall *in vivo* effect may reflect convergence of multiple signalling pathways
 - *Which one is most relevant to use?*

▪ **'Fit for purpose'**

- Suitability for R&D cf. Regulatory environment?
- Practical, meets requirements of clinical sample analysis?

Transferred NAb Assay

Objectives of Transfer Study:

- Confirm suitability of assay for detection of neutralizing antibodies to drug in human serum
- Confirm performance of assays as a pre-requisite for validation
- Demonstrate assay robustness

Target

- Hormone Receptor

Drug type


- Agonist

Mechanism of Action

- Multiple signalling pathways can mediate pleiotropic actions of ligand

Functional Readout

- Cell Proliferation = Convergence of many signalling outputs



Issues with Transfer Assay

Original Assay

- Day to day variability in raw data
- Effected by cell passage / no. days in culture
- Operator Bias
- Required repeat evaluations

Revised Approach

- Analysis of raw data compared with 'normalised' data:
 - ✓ Corrected for day to day variability between runs
 - ✓ Allowed main validation criteria to be met
 - ✗ Unacceptable wide acceptance criteria for positive control

Conclusion

Validation failed due to
inadequate robustness

Agreed Transfer Assay Re-work

- **Modify assay to overcome main obstacles encountered:**

- (1) use of cells in continuous culture
- (2) use of complex downstream readout

- (1) Identified alternate cell line


- Engineered, specific responsive cell line
- No sub-culture requirement; immediate use from liquid N₂

- (2) Functional readout of ‘upstream’ signaling

Direct, Specific, Less Complex

JAK2 = Pivotal Role

Suitable readout for Mechanism of Action



Adapted from Zhu *et al* Cellular Signalling 2001; 13(9); 559

Continuous sub-culture cells

Revive Cell vial



Sub-culture minimum **2 weeks**



Harvest cells; spin & re-suspend x2



Cell count; volume required



Prepare **off-line plate**;
pre-incubate 60min



Incubate Cells + 'reagents'
48 hr (5%CO₂/37°C)



Detection



Read Luminescence

17 days

Immediate use cells

Revive Cell vial



Resuspend, count, **plate**



Incubate **~22 hr** at 5%CO₂/37°C



Prepare **off-line plate** 'reagents'



Incubate Cells + 'reagents'
3 hr (25°C)



Detection



Read Luminescence

2 days

Assays Comparison: Key Parameters

Selected Validation Parameters	Transfer Assay Cell proliferation	Revised Assay 'Upstream' Signalling
Sensitivity to Drug Stimulation	++	+
Minimum Required Dilution	1% serum +	2% serum ++
Screening Assay	Results normalised to ratio	Raw results
Acceptance criteria	Range \pm 30% +	Range \pm 20% ++
Sensitivity to Positive Control	++	+
Inter-assay Precision	+	++
Short-term Stability	Fail?	Pass
Drug tolerance	++	+
Robustness	Fail	Pass

Comparing and Ensuring “Fit for Purpose”

CRITERIA	Transfer Assay: Cell proliferation	Revised Assay Signalling
Cell line / Target Receptor	Mouse Mutated /humanised rat	Human Full length human
Sub-culture required	Minimum 2 week	None
Stability of S/N	Day to day variability	Stable
Impact of sub-culture duration	Variability in sensitivity	Sub-culture not required Uniform cell ‘age’ / condition
Robustness	Inadequate (likely to cause failure of acceptance criteria to be met in SA)	Acceptance criteria meet for sensitivity & system suitability
Intended Use <ul style="list-style-type: none"> ▪ Start up time ▪ Planning Required ▪ Data turnaround time 	Minimum 2 weeks Advanced- adequate flasks for no. plates needed Within 4 weeks	Immediate Stock of vials in liquid N ₂ Within 2 weeks

Summary

- Transfer may commence from different starting points
- Many challenges due to large number of variables involved
- Assay specific pitfalls
- Tried & tested solutions

Critical Factors for Assay Transfer Success

Regular communication & forward planning

Clearly defined objectives & tailored criteria

Continuous monitoring to ensure purpose built approach

Acknowledgements

- David Shaw
- Magda Kieltyka
- Emma Hodgson
- Ray Briggs